

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Escort™ IV Transfection Reagent - Lipid reagent for transient and stable transfection of mammalian and insect cells.

Product Number:

L3287

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Sterility

Conforms

After one week, negative growth on agar plates and in antibiotics free growth medium.

Suitability

Suitable

This lot tested and found to successfully transfect DNA.

Specification: PRD.1.ZQ5.10000044142

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

