

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Monoclonal Anti-γ-Tubulin antibody produced in mouse – clone GTU-88, purified from hybridoma cell culture

Product Number: T5326

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.0 - 1.4 mg/ml
WB-Cell Line/Tissue Extract	HS-68 Human Fibroblast or Chicken Fibroblast
WB Working Concentration	1 - 2 µg/ml
ICC/IF Cells	Hela Cells
ICC/IF Fixation Reagent 1	Methanol
ICC/IF-Fixation Reagent 2	Acetone
ICC/IF-Working Concentration	1 - 2 µg/ml
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000015076

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.