

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**Product Name:
Neutralization Solution BProduct Number: **N3910**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless to Light Yellow
Appearance (Form)	Liquid
No insolubles	
Protein by Biuret	22- 38 mg/ml
Protein by Biuret	2.2- 3.8 %
DNase, Exonuclease Detection	Pass
No More than 90% degradation of the 23,130 bp band from 1 ug of a Lambda DNA Hind III digest marker occurs upon incubation of 25 ul of N3910 at 37 deg. C for one hour	

Remarks:

Specification Date : 08/26/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.