

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Recombinant Trypsin Dimethylated – Proteomics Grade, lyophilized powder, recombinant, expressed in *Pichia pastoris*

Product Number:

EMS0007

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

ug protein/vial

≥ 100

Biuret

units/mg Protein

≥ 10000

Unit Definition: One BAEE unit will produce an A253nm of 0.001 per minute with BAEE as substrate at pH 7.6 at 25 deg C in a reaction volume of 3.2 mL.

Suitability

Suitable

Suitable for use in the digestion of insulin.

Specification: PRD.2.ZQ5.10000087520

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

