

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
SILu™ Lite SigmaMAb Universal Antibody Standard human

Product Number: MSQC4

Storage Temperature: -20 °C

TEST	Specification
mg protein/vial (UV)	≥ 1.0
E1% at 280 nm = 14.0	
Purity by SDS-PAGE	≥ 90 %
Suitability	Suitable
Suitable for use in LC-ESI-MS assays	
Specification: PRD.0.ZQ5.10000053959	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.