

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Invertase from baker's yeast (*S. cerevisiae*) – Grade VII, ≥300 units/mg solid

Product Number:

I4504

CAS Number:

9001-57-4

MDL:

MFCD00131410

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/mg solid

≥ 300

One unit will hydrolyze 1.0 micromole of sucrose to invert sugar per minute at pH 4.5 at 55 deg C.

alpha-Galactosidase

≤ 0.01 %

Recommended Retest Period

1 year

Specification: PRD.0.ZQ5.10000031646

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.