

Product Specification

Product Name:
10x Single Cell Lysis & Fragmentation Buffer

Product Number: L1043

Storage Temperature: -20 °C

TEST	Specification
DNase, Exonuclease Detection	None Detected
NICKase,Endonuclease Detection	None Detected

Specification: PRD.0.ZQ5.10000072628

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.