

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-PUMA/bbc3, N-Terminal antibody produced in rabbit – ~1 mg/mL, affinity isolated antibody, buffered aqueous solution

Product Number:

P4743

MDL:

MFCD04040565

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in phosphate buffered saline, containing 0.02% sodium azide	
Concentration	Conforms
1 mg/mL	
Working Dilution	Conforms
by Direct Immunoblotting	
2-4 ug/ml using human K562 cell lysates. A band of approximately 23 kDa is detected	

Specification: PRD.0.ZQ5.10000087005

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.