

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Amino-terminal FLAG-BAP™ Fusion Protein

Product Number:

P7582

MDL:

MFCD00130921

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
mg protein/ml (Bradford)	1.00 - 10.00
Sensitivity Test	Pass
Minimum 1 ng N-terminal FLAG-BAP detected by dot blot	
Capillary Electrophoresis	≥ 85 %

Specification: PRD.2.ZQ5.10000083466

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

