

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Trypsin inhibitor from Glycine max (soybean) – BioReagent, For use as a marker in SDS-PAGE

**Product Number:** T9767  
**CAS Number:** 9035-81-8  
**MDL:** MFCD00082103

**Storage Temperature:** 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
mg protein/vial (Biuret)	≥ 5
1 vial reconstituted with 1 ml water.	
Suitability	Suitable
Suitable for use as an electrophoresis marker, MW aprox. 20,000 Da	
Recommended Retest Period	Complete
6 years	

Specification: PRD.4.ZQ5.10000007847

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.