

## Product Specification

Product Name:  
β-Glucuronidase Reporter Gene Staining Kit

Product Number: GUSS

Formula Weight: 444,6 g/mol  
Storage Temperature: -20 °C

TEST	Specification
------	---------------

Suitability (Result)	Pass
Suitability (Test)	Enzymatic Assay
Recommended Retest Period	-----
3 Years	

Specification: PRE.1.ZQ5.10000003892

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

