


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polyglycine - mol wt 500-5,000

Product Number: **P8791**
CAS Number: 25718-94-9
MDL: MFCD00166375



Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Solubility (Color)	Light Yellow to Yellow
Solubility (Turbidity)	Clear to Very Slightly Hazy
25 mg/mL, Dichloroacetic Acid	
Water (by Karl Fischer)	≤ 10 %
Degree of Polymerization	9 - 88
by MALDI	
Molecular Weight	500 - 5000
by MALDI	

Specification: PRD.0.ZQ5.10000072197

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.