

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Weigert's iron hematoxylin solution - Part B

Product Number: HT109

MDL: MFCD00284879

TEST	Specification
Appearance (Color) Yellow	Light Yellow to Yellow
Appearance (Form)	Solution
Appearance (Turbidity) Free of particles	Clear
Performance	Conforms
When using Sigma-Aldrich Procedure No. HT30, nuclei stain blue to black. NOTE: Mucin will be stained a deep rose to red.	
Specification: PRD.0.ZQ5.10000031096	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.