

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Acid phosphatase Assay Kit – kit sufficient for 1,000 assays (multiwell plates), kit sufficient for 100 assays (tubes)
Product Number: CS0740

TEST	Specification
Suitability (Result)	Pass
Suitability (Test)	Enzymatic Assay
Suitability (Method)	Technical Bulletin
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000001754

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.