

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Gold chloride solution – 200 mg/dL in deionized water

Product Number:

HT1004



MDL:

MFCD00011322

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

Yellow

Appearance (Form)

Solution

Appearance (Turbidity)

Clear

Essentially free of particles/precipitate

Performance

Conforms

When using with Sigma-Aldrich Procedure No. HT100, fungi stain purplish brown to black.

Specification: PRD.4.ZQ5.10000028001

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.