

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Adenosine 5'-triphosphate disodium salt solution - Crystalline ATP, HPLC purified, aqueous solution for RNA transcription

Product Number:

A6559

MDL:

MFCD00080339



Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Appearance (Turbidity)

Clear

DNase, Exonuclease Detection

None Detected

RNase Detection

None Detected

Purity (HPLC)

≥ 95 %

Concentration

95 - 105 mM

Based on EMM at 259 nm = 15.4

Specification: PRD.1.ZQ5.10000023558

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.