

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Decellularized ECM bioink precursor – from porcine skin, suitable for 3D bioprinting applications

Product Number:

906867

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Color)

White to Beige

ISpec - to be done by sampling for verification purposes

Appearance (Form)

Sponge

ISpec - to be done by sampling for verification purposes

Appearance (Color)

White to Beige

Appearance (Form)

Sponge

(dry)

Product of Supplier

Confirmed

T&R Biofab

Specification: PRD.0.ZQ5.10000103793

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

