

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Laminin from human fibroblasts – cell culture derived, liquid, sterile-filtered

Product Number: L4544

Storage Temperature: -70 °C

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Solution
Sterility	Sterile
SDS Electrophoresis	Conforms
High molecular weight laminin bands detected	
Country of Origin	United States of America
Country of Manufacture	United States of America
Tested For Infectious Agents	Tested
Protein Content	0.3 - 0.6 mg/ml
in TBS, pH 7.4	

Specification: PRD.1.ZQ5.10000013452

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.