

## Product Specification

**Product Name :** Glutathione Peroxidase Assay Kit Sufficient for 100 colorimetric tests  
**Product Number:** MAK437  
**Storage Temperature :** -20°C

Test	Specification
Suitability The kit is suitable for the quantitative determination of glutathione peroxidase activity and evaluation of drug effects on GPX activity in biological samples.	Suitable
Slope	$\geq 0.2$
Assay Linearity R2	$\geq 0.98$

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1203841

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact [mlsbranding@sial.com](mailto:mlsbranding@sial.com)

