

Product Specification

Product Name:

Trifluoroacetic acid - for protein sequence analysis

Product Number: 1.08178

CAS Number: 76-05-1

Formula: C₂HF₃O₂

Formula Weight: 114,02 g/mol

TEST	Specification
Assay (Acidimetric)	≥ 99.7 %
Density (D 20/20); (About 1.53)	
Water	≤ 0.02 %
Infrared Spectrum	Passes Test
Miscellaneous Assay (Boiling Point; About 72°C)	Passes Test

Specification: PRE.0.ZQ5.10000025881

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

