

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

123 bp DNA Ladder – for DNA electrophoresis

Product Number:

D5042

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Suitability

Confirmed

Suitable for use as a DNA marker.

Concentration

Confirmed

0.8 - 1.2 mg/ml

Specification: PRD.1.ZQ5.10000022021

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.