

## Product Specification

**Product Name :** Thrombin from bovine plasma lyophilized powder, 600-2,000 NIH units/mg protein (biuret)  
**Product Number:** T6634  
**CAS Number:** 9002-04-4  
**MDL Number :** MFCD00166053  
**Storage Temperature :** -20°C

Test	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
NIH units/mg protein (Biuret) Unit Definition: Activity is expressed in NIH units obtained by direct comparison to a NIH Thrombin reference standard, Lot K	600 - 2000
Recommended Retest Period 5 years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1152111

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact [mlsbranding@sial.com](mailto:mlsbranding@sial.com)

