

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-phospho-c-Abl (pTyr⁴¹²) antibody produced in rabbit - affinity isolated antibody, buffered aqueous solution

Product Number:

C5240

MDL:

MFCD06798483

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Formulation	Conforms
Solution in Dulbecco's phosphate buffered saline, pH 7.3, 50% glycerol, with 1 mg/mL BSA and 0.05% sodium azide	
Working Dilution	Conforms
by Western Blot Recommended starting dilution 1:1000. The optimal antibody concentration should be determined by user for application.	
The optimal antibody concentration should be determined by user for application.	
Volume	Conforms
100 uL (sufficient for 10 tests)	

Specification: PRD.2.ZQ5.10000074232

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

