

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


Poly(ethyleneimine) solution - average  $M_w$  1800, 50 wt. % in  $H_2O$ 

Product Number:

408700

MDL:

MFCD00084427



TEST	Specification
------	---------------

Appearance (Color)

Colorless to Faint Yellow

May Darken to Faint Yellow Over Time

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Viscosity

320 - 350 m.Pa/s

Assay

48.0 - 52.0 %

Dry Content

Molecular Weight

1500 - 2000

by GPC

Specification: PRD.4.ZQ5.10000029373

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

