

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Anti-MAP2 antibody produced in rabbit – ~1 mg/mL, affinity isolated antibody, buffered aqueous solution

Product Number: **M3696**

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Solution
in PBS containing 0.02% sodium azide	
Concentration	1 - 1 mg/ml
Western Blot	Dilution Range 1:500 - 1:1000
Immunohistochemistry	Dilution Range 1:100 - 1:200

Specification: PRD.0.ZQ5.10000020548

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.