

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Iodomethane - contains copper as stabilizer, ReagentPlus®, 99%

Product Number:

I8507

CAS Number:

74-88-4

MDL:

MFCD00001073

Formula:

CH₃I

Formula Weight:

141.94 g/mol

Storage Temperature:

2 - 8 °C



TEST

Specification

Appearance (Color)

Colorless to Faint Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.5 %

Color Test

≤ 100 APHA

Note

Confirmed

Contains Copper as Stabilizer

Recommended Retest Period

1 year

Specification: PRD.4.ZQ5.10000005109

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

