

Product Specification

Product Name : Lipid Mixture 1, Chemically Defined liquid, sterile-filtered, BioReagent, suitable for cell culture
Product Number: L0288
Storage Temperature : +2°C to +8°C

Test	Specification
Appearance (Turbidity)	Clear to Hazy
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Endotoxin Level	≤ 10EU/ml
Sterility by USP guidelines	Negative
Cell Culture Test at 10ml/L	Pass
Cell Line	Attached Cell Line
Cell Line	Attached Cell Line
Cell Line	Attached Cell Line

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1206027

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

