

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


OptiPrep™ Density Gradient Medium - used for cell and subcellular organelle isolation

Product Number:

D1556

MDL:

MFCD00867965



TEST

Specification

Bulk Density	1.319 - 1.321 g/ml
Endotoxin Level IU/mL	≤ 1
Sterility	Pass
Trademark	Confirmed
OptiPrep is a trademark of Serumwerk Bernburg AG	

Specification: PRD.3.ZQ5.10000033831

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

