

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Latex beads, deep blue dyed - 0.24 µm mean particle size, aqueous suspension, solids 10 %

Product Number:

L1273

MDL:

MFCD00131492

TEST

Specification

Appearance (Color)

Light Blue to Very Dark Blue

Appearance (Form)

Suspension

Particle Diameter

0.23 - 0.25 micron

Standard Deviation

≤ 0.02 micron

Specification: PRD.1.ZQ5.10000025921

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

