

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
LB Agar Plates – pre-poured agar plates for E. coli

Product Number: L5542

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance	Conforms
Amber agar plates	
Sterility	Pass
Performance Testing	Pass
Growth promotion and Antibiotic or Amino acid sensitivity	

Specification: PRD.0.ZQ5.10000037317

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.