

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(Bisphenol A-co-epichlorohydrin), glycidyl end-capped - average $M_n \sim 377$

Product Number:

405493

CAS Number:

25036-25-3

MDL:

MFCD00197891



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Gardner Color	≤ 1
Epoxide Equivalent Weight	185.0 - 192.0 g/eq
Viscosity	100 - 150
At 25°C, Poise	

Specification: PRD.0.ZQ5.10000107182

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

