

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
ProtectRNA™ RNase Inhibitor 500× Concentrate – RNase inhibitor for in situ hybridization assays

Product Number: R7397

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Red to Very Dark Red
Protect RNA Dilution	1:500
Suitability	Pass
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000017574

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.