

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
TissueFab® bioink Conductive - Vis/405 nm

Product Number: 915963

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Dark Grey to Black
Appearance (Form)	Viscous Liquid or Gel
pH	6.5 - 7.5
Bioburden	≤ 5 CFU/g
(Fungal)	
Bioburden (Total Aerobic)	≤ 5 CFU/g
Trademark	Note
TissueFab is a registered trademark of Sigma-Aldrich Co. LLC	

Specification: PRD.0.ZQ5.10000107137

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

