

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hyaluronic Acid Binding Protein bovine – lyophilized powder

Product Number: H0161

MDL: MFCD02099073

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
SDS Electrophoresis	Consistent with Past Lots
Band Pattern by SDS Electrophoresis	
Purity by SDS Electrophoresis	Conforms
> 90%	

Specification: PRD.0.ZQ5.10000014887

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.