

Product Specification

Product Name : p-Acetamidophenyl β-D-glucuronide sodium salt >98% (TLC)
Product Number: A4438
CAS Number: 120595-80-4
MDL Number : MFCD00153898
Storage Temperature : -20°C
Molecular Formula: C₁₄H₁₆NNaO₈
Formula Weight: 349.27 g/mol

Test	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity) 25 mg/mL, H ₂ O	Clear
Water (by Karl Fischer)	≤ 6%
Sodium (anhydrous)	6.3 - 7.0 %
Purity (TLC)	≥ 98%
IR Spectrum	Conforms to Structure
Recommended Retest Period 4 years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1158932

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

