

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

PCR 50 – 2,000 bp Marker – for electrophoresis of PCR fragments

Product Number:

P9577

Storage Temperature:

-20 °C

TEST

Specification

Agarose Electrophoresis

Conforms

Resolved 8 bands consistent with 50, 150, 300, 500, 750, 1000, 1500, and 2000 base pairs.

Specification: PRD.1.ZQ5.10000013360

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.