

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Gold nanowires – diam. \times L 30 nm \times 4,500 nm, dispersion (H₂O), contains CTAB as stabilizer**Product Number: 716944**

Formula Weight: 196.97 g/mol

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Colorless to Light Orange
Appearance (Form)	Conforms
Liquid with White to Off-White Flocculant (Precipitate)	
Size	24 - 36 nm
Axial Diameter	
Size	3600 - 5400 nm
Longitudinal Size	
Concentration	\geq 50 µg/ml
Recommended Retest Period	-----
1 Year	

Specification: PRD.0.ZQ5.10000021798

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.