

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-Albumin antibody produced in rabbit – IgG fraction of antiserum

Product Number: A0433

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	10.0 - 20.0 mg/ml
Antibody Concentration(Method)	Quantitative Precipitin Assay
Antibody Concentration	4.0 - 6.0 mg/ml
ELISA - Titer	1:25,000-1:50,000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000015672

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.