

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

GenElute™ Bacterial Genomic DNA Kits – sufficient for 350 purifications

Product Number:

NA2120

TEST	Specification
------	---------------

Suitability

Suitable

The average yield must be at least 70% of the average yield of the prequalified control solution.

Level of degradation seen pBR322 DNA isolated with test solution and digested with HindIII should be comparable to the control when examined by agarose gel electrophoresis.

Specification: PRD.0.ZQ5.10000067088

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.