

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Black phosphorus - ≥ 20 micro gram/ml in NMP

Product Number:

902896



TEST	Specification
------	---------------

Appearance (Color)

Yellow to Very Dark Yellow

Appearance (Form)

Liquid

Infrared Spectrum

Conforms to Structure

Phosphorus (P)

 $\geq 20 \mu\text{g/ml}$

Recommended Retest Period

1 Year

Specification: PRD.2.ZQ5.10000099037

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

