

Product Specification

Product Name : Staphylococcal enterotoxin A from Staphylococcus aureus
Product Number: S9399
CAS Number: 642595-84-4
Storage Temperature : Cooler/Refrigerated

Test	Specification
Appearance (Form)	Lyophilized Powder
Appearance (Color)	White to White w/ Yellow Cast
Solubility (Solvent)	Water
Solubility (Conc)	1.00 - 2.00 mg/ml
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless
Protein Content (Method)	Lowry
Gel Electroph. (Gel)	SDS-Page
Gel Electroph. (Stain)	Coomassie
Gel Elec. (Compar. to Cont)	Similar
Gel Electroph. (Description)	1 Major Band
Recommended Retest Period 3 Years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1098157

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

