

Product Specification

Product Name :

EX-CELL® CD HEK293 Viral Vector Medium Chemically defined, animal component-free, without L-glutamine, liquid, sterile-filtered, suitable for cell culture

Product Number: 14385C

Storage Temperature : +2°C to +8°C

Test	Specification
Appearance (Turbidity)	Clear
Appearance (Form)	Solution
pH	7.0 - 7.2
Osmolality (mOsm/kg)	290 - 310
Sterility by USP Guidelines	Conforms
Endotoxin Level Record	
Site of Manufacturing	Broadway

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1266502

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

