

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

CleanAmp™ dNTP – Hot-start dNTP mix for improved PCR

Product Number:

DNTPCA1

Storage Temperature:

-20 °C

TEST	Specification
Kit Comp (Component 1)	-----
CleanAmp(TM) dATP	
AX-HPLC	Pass
RP-HPLC	Pass
31P-NMR	Pass
1H-NMR	Pass
Kit Comp (Component 2)	-----
CleanAmp(TM) dGTP	
AX-HPLC	Pass
RP-HPLC	Pass
31P-NMR	Pass
1H-NMR	Pass
Kit Comp (Component 3)	-----
CleanAmp(TM) dCTP	
AX-HPLC	Pass
RP-HPLC	Pass
31P-NMR	Pass
1H-NMR	Pass
Kit Comp (Component 4)	-----
CleanAmp(TM) dTTP	
AX-HPLC	Pass
RP-HPLC	Pass
31P-NMR	Pass
1H-NMR	Pass
Functional Assay	Pass
CleanAmp(TM) dNTP Mix	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

CleanAmp™ dNTP – Hot-start dNTP mix for improved PCR

Product Number:

DNTPCA1

Storage Temperature:

-20 °C

TEST

Specification

Specification: PRD.0.ZQ5.10000063845

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.